

# CRNA TRAINING KITS

LATEX FREE

Ideal for anesthesia professionals in training, CRNA training kits contain essential tools for residency work and beyond. Each kit contains a peripheral nerve stimulator, monoscope, four piece chest set and durable canvas pouch.



CRNA Professional Training Kit



CRNA Student Training Kit



CRNA Basic Training Kit

## CRNA TRAINING KITS

- Essential diagnostic tools for anesthesia students
- Three kit options feature either a SunStim™, SunStim™ Plus or SunStim™ Pro nerve stimulator
- Priced right for student budget
- One convenient purchase rather than separate product orders
- Durable canvas pouch keeps all tools close at hand

ITEM	DESCRIPTION	PK
6-1030-53	CRNA Professional Training Kit	1
6-1030-52	CRNA Student Training Kit	1
6-1030-50	CRNA Basic Training Kit	1

## CRNA TRAINING KITS

### COMPONENTS:



Basic Kit



Student Kit



Professional Kit

## PERIPHERAL NERVE STIMULATOR

### SUNSTIM™

- Twitch: 1 pulse per second
- Tetanus: 100 Hz (100 pulses/s)
- Built-in train of four (TOF): 4 pulses in a period of 2 seconds
- Stimulator output current is fully adjustable from 0-70 mA

### SUNSTIM™ PLUS, SUNSTIM™ PRO

- Double burst (DBS): Two 60-ms bursts of 50 Hz separated by 0.75s
- Twitch: 1 pulse per second
- Tetanus: 50 Hz (50 pulses/s) or 100 Hz (100 pulses/s)
- Built-in train of four (TOF): 4 pulses in a period of 2 seconds
- Stimulator output current is fully adjustable from 0-70 mA
- Stimulus pulse: square wave monophasic pulses (200- $\mu$ s duration)



## MONOSCOPE

- Designed specifically to address sound transmission challenges
- Once compressed foam earpiece is inserted, it expands slowly for a perfect, low pressure seal with the canal wall
- Sealing action enables the accurate transmission of low volume sounds
- 4 ft Tubing



## CHEST PIECES

- Durable stainless steel construction
- Pediatric chest pieces are hollow and lightweight, with rounded, smooth chest surface, ideal for use with neonates
- Solid Wenger chest pieces feature a relatively deep sound chamber
- Luer slip male connector
- Comes with 8 double adhesive disks